


FEATURES	ZOLL AED PLUS	HEARTSINE SAMARITAN PAD 300 P	HEARTSTART FIRST AID DEFIBRILLATOR
Size (H x W x D)	13 cm x 24 cm x 29 cm	5 cm x 20 cm x 18 cm	7 cm x 21 cm x 19 cm
Weight	3.1 kg	1.1 kg	1.5 kg
Battery capacity (Standby)	Up to 5 years (typical)	Up to 3.5 years (typical)	Up to 4 years (typical)
Battery capacity (Shocks - use)	300 Shocks (13 hours)	30 shocks (6 hours)	90 shocks (3 hours)
Battery type	Type 123A Photo Flash	18 V LiMnO2 disposable	9 V DC LiMnO2 disposable
Pads shelf life	Up to 5 years (typical)	Up to 3.5 years (typical)	Up to 2 years (typical)
Use in children (< 25 Kg)	YES (Pedi Padz)	YES (Pediatric-Pak)	YES (Smart Pads cartridge)
CPR feedback	YES (Adults)	NO	NO
CPR Voice Coaching	YES	YES	YES
Water / Dust protection (IP)	55 (Dust protected / Water jets)	56 (Dust protected / heavy seas)	21 (> 12 mm / dripping water)



be prepared...
now you can help save a life



We cannot rely on an Ambulance or a Rescue Helicopter being only minutes away. These resources can take time to arrive - in the city, in the wilderness or on the water.

Seconds count! The only proven effective treatment for ventricular fibrillation (VF) is quick defibrillation. When someone has a Sudden Cardiac Arrest (SCA), survival may be significantly increased if early defibrillation is administered within the first few minutes. The graph to the left shows the relationship between sudden cardiac arrest survival rate & time to defibrillation*.

*Reference: AHA 2000 Guidelines

You never know when sudden cardiac arrest will happen. Any time, any place. At home, in the office, at work, the club or whilst travelling. **Seconds count!**

CONTACT US:

- Telephone: 1300 729 575
- Facsimile: 02 8003 9875 ■ Email: info@defibshop.com.au
- Defib Shop PO Box 367, Willoughby NSW 2068



www.defibshop.com.au

www.defibshop.com.au





HAVE THE PEACE OF MIND KNOWING THAT YOU CAN SAVE A LIFE EVEN BEFORE HELP ARRIVES. EASY TO USE BY ANYONE, ANY TIME, ANYWHERE!

- Seconds count! Every minute that passes without a defibrillator decreases your chance of survival
- An Ambulance or Rescue Helicopters may be delayed in arriving to help
- Defibrillators can increase your chance of survival in cardiac arrest from approx. 10% to 60%
- Easy to learn and use
- Visual and/or audio prompts to assist rescuer
- Compact, lightweight, portable and easy to transport
- Compliant with Australian Resuscitation Council guidelines

ZOLL AED Plus Real CPR Help®

The ZOLL AED Plus has proven popular with schools, clubs, recreation and holiday parks and workplaces.

- Low ongoing costs
- Real CPR Help® provides feedback on quality of CPR
- Uses readily available 'camera' batteries providing up to 300 shocks
- Visual and verbal prompts to assist the rescuer
- Single piece defib pads to make application easy
- Suitable for use on paediatric patients with optional defib pads.



HeartStart First Aid Defibrillator

The HeartStart First Aid Defibrillator is popular in public places, working environments and clean and dry workplaces.

- CPR coaching to support the rescuer
- Lightweight and compact
- Paediatric cartridge optional
- Metronome to assist maintain CPR rate.



Heartsine Samaritan 300P

The Samaritan 300P is popular where weight and size matter! For example hiking, camping and bushwalking, aviation and boating. Recommended when you must carry a defibrillator with you.

- Portable and lightweight at just 1.1kg including batteries and defib pads
- Integrated battery and defib pad pack
- CPR coaching to support the rescuer
- IP rating of 56 - protected for wet and dusty environment
- Paediatric Smart Pads cartridge optional.

