

## HeartSine Products and Accessories

July 2009

### Products: HeartSine samaritan PAD and Trainer



#### HeartSine samaritan® PAD systems (model 300P)

Includes 7 year warranty, IP56 Certification, user manual, soft carry-case, quick reference card. Also includes a Pad-Pak™ (see list below for available configurations). Pad-Pak™ includes battery and 1 pair of electrode pads, and carries 3.5 year expiration.

**HeartSine samaritan PAD with 1 Adult Pad-Pak™**, Part# PAD-BAS-US-05

**HeartSine samaritan PAD with 2 Adult Pad-Paks™**, Part# PAD-SYS-US-05

**HeartSine samaritan PAD with 1 Adult Pad-Pak™ & 1 Pediatric-Pak™**

Part# PAD-PED-US-05




#### HeartSine samaritan® PAD TRAINER







Includes: 2 year warranty (unit & charger), Trainer-Pak™ (1 yr. warranty), user manual, soft carry case, 5 sets Training pads. Also includes: Pause button, 8 training scenarios (4-AHA/ERC 2005 and 4-2000 guidelines), 120-240 Volt system with rechargeable battery.

Part# TRN-SYS-US-05

### Accessories: HeartSine samaritan PAD (models 300/300P)

See list below for replacement accessories available for the HeartSine® samaritan® PAD. Pad-Pak™ expiration dates are clearly marked underneath each unit. Replace before expiry date!

 <p><b>Adult Pad-Pak™ (gray)</b> Part# PAD-PAK-01</p>  <p><b>Pediatric-Pak™ (pink)</b> For children &lt; 8 years of age or &lt; 55 lbs Part# PAD-PAK-02</p>	 <p><b>Wall Cabinet, Custom</b> Alarmed, dimensions 12"H x 11"W x 5" D, weight 5 lbs. Shown with custom glass etching. Part# PAD-CAB-04</p>
 <p><b>Prep Kit</b> Includes mask, razor, scissors, alcohol wipe Part# PAD-ACC-17</p>	 <p><b>Wall Cabinet</b> Alarmed, with custom glass front, dimensions 13"H x 13.5"W x 7"D Part# PAD-CAB-03</p>
 <p><b>SAVER® EVO Software and Cable Kit</b> Event recovery or print software for PAD and AED; with ECG, shock information, CPR intervals, etc. Also use to switch PAD/AED language. Includes USB cable. Requires Windows 2003 or XP, &amp; USB port. Part# PAD-ACC-01</p>	 <p><b>Wall Cabinet</b> Alarmed, with standard glass front, dimensions 13"H x 3.5"W x 7"D Part# PAD-CAB-01</p>

 <p><b>Universal Software Upgrade CD and Cable Kit, for upgrading models 300/300P</b> Software also available via heartsine.com Part# PAD-ACC-04-03</p>	 <p><b>PAD Wall Bracket</b> Premium Wire wall bracket for convenient storage of your PAD in an office or vehicle Part# PAD-CAB-02</p>
 <p><b>USB Cable</b> Use with Saver EVO data management software and Universal Upgrades. Software available via heartsine.com or in Kits. Part# PAD-ACC-02</p>	 <p><b>PAD Hard Case</b> Watertight carry case for further protection against drops and for use in a vehicle, swimming pool, marine and camp applications, etc. Case size = 17" x 14" x 7". Weighs 9 lbs. Part# PAD-BAG-02</p>
 <p><b>AED Wall Sign – Triangular</b> Part# PAD-ACC-05</p>	 <p><b>Replacement Soft Carry-Case</b> (AED not included), Part# PAD-BAG-01</p>

## Accessories: HeartSine samaritan PAD Trainer



### Trainer Defibrillation Pads

Part# TRN-ACC-02 (pkg/10)      Part# TRN-ACC-03 (pkg/25)



### Replacement Battery Charger

Part# TRN-ACC-01



### Replacement Trainer-Pak

Part# TRN-PAK-01

## Accessories: HeartSine samaritan AED (models sam001/002/003)



### SDRU (software data recovery unit)

Part# SDU-003



### Data-Pak Battery (32 MB)

Part# SBP-001



### AED Electrode (Pair)

Part# SDE-201 (pkg/1); SDE-210 (pkg/10)

### Corporate Headquarters:

HeartSine Technologies, Inc.  
121 Friends Lane, Suite 400  
Newtown, PA 18940  
Tel: (215) 860 8100  
Fax: (215) 860 8192  
Web: heartsine.com

### Manufacturing Facilities

HeartSine Technologies, Ltd  
Canberra House, 203 Airport Road West  
Belfast, Northern Ireland BT3 9ED  
Tel: +44 (0) 28 9093 9400  
Fax: +44 (0) 28 9093 9401

### Distributed by:



The products described in this brochure all meet the applicable European Medical Directive requirements.



**HeartSine®**

*Lifesaving Technology  
Designed for the Real World.  
Since 1998.*